

Beer calibration sample set / Bierkalibrierproben Satz



	Unit	Sample ID / Probennummer										Range / Bereich		Unit
	Einheit	BCS0001	BCS0002	BCS0003	BCS0004	BCS0005	BCS0006	BCS0007	BCS0008	BCS0009	BCS0010	Min	Max	Einheit
pH value (EBC) pH-Wert (EBC)		4,73	4,22	4,45	4,60	4,12	4,09	4,34	4,42	4,27	3,86	3,86	4,73	
Color (EBC) Farbe (EBC)	EBC	11,90	16,57	105,30	8,80	106,33	7,13	30,80	12,40	7,47	36,25	7,13	106,33	EBC
Bitter unit Bittereinheiten (EBC)	IBU	19,63	13,25	20,52	20,70	26,58	14,35	21,00	21,67	20,50	14,12	13,25	26,58	IBU
Density Beer Dichte Bier	g/cm³	1,00779	1,00646	1,00780	1,02499	1,00583	1,00356	1,00710	1,00853	1,00644	1,00749	1,00356	1,02499	g/cm³
Apparent extract Scheinbarer Extrakt	% w/w	2,46	2,12	2,47	6,78	1,96	1,38	2,28	2,65	2,12	2,38	1,38	6,78	% w/w
Real extract Wirklicher Extrakt	% w/w	5,10	5,49	4,12	6,87	3,49	3,00	4,24	4,22	3,79	4,03	3,00	6,87	% w/w
Alcohol by weight Alkohol	% w/w	6,03	7,99	3,47	0,31	3,29	3,65	4,34	3,76	3,79	3,68	0,31	7,99	% w/w
Alcohol by volume Alkohol	% v/v	7,70	10,19	4,77	0,41	4,20	4,64	5,54	4,81	4,84	4,70	0,41	10,19	% v/v
Original gravity Stammwürzegehalt	% w/w	16,50	20,27	11,39	7,49	9,95	10,15	12,63	11,53	11,18	11,20	7,49	20,27	% w/w
Apparent attenuation Vergärungsgrad scheinbar	% w/w	85,10	89,50	78,30	9,44	80,30	86,40	81,91	77,00	81,00	78,70	9,44	89,50	% w/w
Real attenuation Vergärungsgrad wirklich	% w/w	70,90	75,10	65,20	8,59	66,10	71,50	67,92	64,80	67,40	65,30	8,59	75,10	% w/w

For further info please check the certificate of analysis / Weitere Informationen entnehmen Sie bitte dem Analysenzertifikat